

**OsteoSys**

better products for better life

New portable QUS for everyone

# Beetle



3.8kg



5Sec.



NFC



GPS



Bluetooth

- Portable wireless QUS(Quantitative Ultrasound) bone mineral densitometer
- Compact size, utmost lightness and enough space for calcaneus measurement
- Connected to mobile phone by bluetooth
- Easiest & quickest measurement with one-touch
- Efficient result reporting and data exporting by special app.
- Location tracking by GPS
- Multi battery system : Built-in battery interchangeable with backup battery( $\geq 5V$ , 2A)
- Temperature compensation function
- LED and audio guidance



New QUS for Everyone



BeeTLe is the one and only QUS bone mineral densitometer with NFC and GPS functions.  
The measurement can be briefly done by one-touch, and the whole measuring process is operated by BeeTLe's special app. on your mobile phone.



3.8kg



5Sec.



NFC



GPS



Bluetooth

- QUS (Quantitative Ultrasound)
- Measurement site : Calcaneus
- Measurement time : 5sec.
- Reproducibility : BUA -  $\leq 1.5\%$  C.V.  
SOS -  $\leq 0.2\%$  C.V.  
BQI -  $\leq 1.5\%$  C.V.
- Measured parameter : T-score, Z-score, BQI, BUA, SOS,  
Ratio of T-score and Z-score
- Dimension : 249 × 337 × 146mm
- Weight : 3.8kg
- Pediatrics / Trend report / DICOM & PACS
- NFC / GPS / Bluetooth / App. / LED & Sound



**OsteoSys**

OsteoSys Co., Ltd.

9F, 903 JnK Digital Tower, 111 Digital-ro 26, Guro-gu, Seoul, Korea Tel. +82.2.6124.5900 Fax. +82.2.6124.5958 [www.osteosys.com](http://www.osteosys.com)

